


Note: Electrical requirements are 8 Amps / 115 Volts/ 60Hz/ 3 Pole GFI Receptacle.



All dimension figures
in parentheses are
in millimeters.



SLANT BACK STATIONARY WHIRLPOOL SB-100-S

NOTES:

- | | |
|-------------------------|---------------------------------------|
| 1. Slant Back Whirlpool | 5. Turbine Raising & Lowering Device. |
| 2. Filler Spout | 6. Drain & Overflow Assembly |
| 3. Thermometer. | 7. Floor Mounting Bracket Assembly |
| 4. 1/2 HP Pump Motor. | |