

**6 Drawer MRI-Safe Emergency Cart Specialty Package**  
**Item# GME-B6RMHRF**

The emergency cart is designed and built for safe use in the strong magnetic fields present in magnetic resonance imaging (MRI) suites.

**Features:**

- IV Pole
- Cardiac board
- O2 holder
- Utility Hooks (4)
- Full drawer trays with adjustable dividers (2)
- Tall cabinet with breakaway locking system
- Vinyl-padded full wrap-around steel bumper
- 5" MR-safe single-wheel caster set, full swivel – two with brake
- Stainless steel pull-out shelf
- Pre-threaded accessory mounting holes
- Low ferrous ball-bearing drawer slides
- Replaceable plastic top with molded-in push handles
- Cart is delivered fully assembled
- Latex free
- Certified MR-safe by independent testing authority (copies of report available)
- MR-safe accessories available
- Hard-baked powder coat finish
- Shelf painted gloss white with MRI sticker



**Specifications:**

- Height: 45.5"
- Width: 32"
- Depth: 25"
- IV pole extends to 68.75in H
- Weight: 144 lbs



**Drawer Configuration:**

- (1) – 4in
- (2) – 3.25in
- (2) – 6.5in
- (1) – 9.75in

Internal drawer dimensions: 23.9in W x 17in D

**Manufacturer warranty: 12 yr**



Toll Free 1.877.706.4480 | Local/ International 713.706.4480 | Fax 713.706.4161  
[www.MSEC.co](http://www.MSEC.co)