

## Lite Mobile Emergency Crash Cart with Accessories Item# MS3C-110CR

This treatment cart offers large capacity and flexible organization of treatment supplies. This cart includes an expansive work area and an additional slide out work surface for additional workspace. Durable caster position ensure stability while making the cart easy to move and control.

**Quick Ship Program, ships in 2 days! Free Dock to Dock shipping, Continental US!**

### Features:

- Made of ABS plastic, easy to clean and corrosion-proof
- All-in-one table top with raised sides to eliminate items rolling off
- Built-in handrail for easy maneuvering
- ABS table top with transparent dust-proof pad
- Sliding side shelf for writing
- 5 ABS drawers slide smooth and silent with centralized lock and plug-in labels to avoid dust
- (3)Flexible material storage bins that can be used for gloves, etc.
- (4) Anti-Static castors, two locking, designed to eliminate debris from wrapping the caster
- Drawer dividers for arbitrary partitioning
- (2)Trash bins
- Telescoping IV pole
- Defibrillator Shelf
- C.P.R board
- Oxygen tank holder



### Drawer Configuration:

- (2) Small drawers: 16.5in (W) x 14in (D) x 3.5in (H)
- (2) Medium drawers: 16.5in (W) x 14in (D) x 5in (H)
- Large drawer: 16.5in (W) x 14in (D) x 8in (H)

### Dimensions:

- 19 in (D) x30 (W) in x 38in (H)

**Manufacturer warranty: 1 year**



Toll Free 1.877.706.4480 | Local/ International 713.706.4480 | Fax 713.706.4161  
[www.MSEC.co](http://www.MSEC.co)