



toll free: 1-877-706-4480  
international/local: 713-706-4480  
fax: 713-706-4161  
[www.MS3C.co](http://www.MS3C.co)

## Lite Mobile Medication Cart w/ Accessories

Item# MS3C-200M

This medication treatment cart offers large capacity and flexible organization of medications. This cart includes an expansive work area. Durable caster positioning ensures stability, while making the cart easy to move and control.

### Features:

- Made of ABS plastic and steel-plastic, easy to clean and corrosion-proof
- ABS table top with transparent dust-proof liner mat
- (2) Key locks for safe storage, one lock for each side
- Steel guardrail stabilizes items during transport
- Sliding side shelf for writing
- (6) ABS drawers on each side (12 drawers total) slide smooth and silent with centralized lock and plug-in labels to avoid dust
- (4) 5" Heavy Duty Anti-Static castors, two locking, designed to eliminate debris from wrapping the caster
- Two keys included
- (1) Chart Storage Bin
- (2) Trash bins



### Drawer Configuration:

- (10) Small drawers: 22in (W) x 9in (D) x 3in (H)
- Large drawer: 22in (W) x 9in (D) x 15in (H)

### Dimensions:

- 34 in (L) x 22.5 (W) in x 42in (H)
- Shipping Dims: 37 in (L) x 25 (W) in x 44in (H)
- Approximate weight: 85 lbs.

**Manufacturer warranty: 1 year**