

Low Height Bariatric Treatment Table

Item# 6916CLT

This low-height bariatric table features contemporary design, convenient storage and easy clean laminate shelf. It is available in two base colors and nine vinyl upholstery colors to enhance and complement any decor.

Features:

- Lower height makes patient access easier
- Plywood top and solid hardwood legs
- 6-leg design with triple bolted corner legs
- Built-in middle leg levelers
- Steel reinforced frame
- Laminate full length shelf with H-braces
- Pillow included
- Seamless rounded corners
- Two-inch firm foam padding
- Pre-threaded inserts make replacing a damaged top quick and easy
- Easy-clean and sanitize all laminate shelf
- Two colors of hardwood legs
- 9 colors of premium, stain-resistant, woven, knit-backed vinyl



STANDARD VINYL UPHOLSTERY—Corresponding color code must accompany color selection when placing an order.



Specifications:

- Dimension: 78in L x 26in H x 30in W
- Load capacity: 600lbs

Manufacturer warranty: 5 years



Toll Free 1.877.706.4480 | Local/ International 713.706.4480 | Fax 713.706.4161
www.MSEC.co