



toll free: 1-877-706-4480
international/local: 713-706-4480
fax: 713-706-4161
www.MS3C.co

EMS Evacuation Stair Chair **Item# MS3C-300TS**

This durable and easy to operate Emergency Evacuation Stair Chair is designed to aid in stairway descent during emergency situations or in case of power failure. Designed for an almost instantaneous deployment and operation on any surface, this evac stair chair is a must safety device for multi story residential or commercial buildings.

Ships Free, Continental US Only.

Features:

- Dual, smooth rubberized tracks to aid care providers with controlled descent on any surface; carpet, wood, metal, and concrete.
- Back of the evacuation chair is designed with two folded handles that lock into position allowing operators multiple grip options for safe handling and easy steering
- Two flexible uplift rods under the front of the stretcher to aid in lifting
- Made of high-strength aluminum alloy materials
- Durable seat material is easy to clean
- Easily folds into a compact design making it easier for storage and transport
- Light-weight, easily-carried
- Safety shoulder and lap belts Included

Optional Accessories:

- Dust Cover
- Wall mount brackets

Specifications:

- Working load capacity: up to 400 lbs.
- Upright dimensions: 51.5in(L) x 19.5in (W) x 61in (H)
- Packing size: 42in x 21in x 10in
- Net weight: 31 lbs

