



toll free: 1-877-706-4480  
international/local: 713-706-4480  
fax: 713-706-4161  
[www.MS3C.co](http://www.MS3C.co)

## EMS Series Ambulance Stretcher

### Item# MS3C-1000

The MS3C-1000 EMS stretcher is an affordable, lightweight, durable stretcher that's easy to operate. The stretcher can be folded into a transport chair permitting access to all those confined spaces such as corridors with narrow angles or elevators.

#### Features:

- Ideal solution for patients who cannot be subjected to long waiting time in a flat position
- Both ends of the stretcher can be adjusted to 12 various positions
- Mattress is made of 2in waterproof foamed cushion
- Made of high-strength, lightweight aluminum alloy material
- Height of the backrest and knee rest is adjusted with a pneumatic spring
- Only one person required to operate the stretcher when loading it into an emergency vehicle
- Mechanical folding legs are controlled by the handles on both sides and automatically lock in place when the undercarriage is fully unfolded
- Stretcher comes standard with vest-style safety belts
- Easy to disinfect and clean

#### Specifications:

- Upright dimensions: 77in x 22in x 33.5in
- Collapsed dimensions: 76.8in x 22in x 9.8in
- Full upright chair dimensions: 54in H x 46in L x 22in W
- Maximum Vehicle Loading Height: 26.5 in
- Shipping size: 78in x 24.4in x 10.6in
- Load capacity: 350 lbs

**Manufacturer warranty: 1 year**

