

# MS3C-100S Lightweight Automatic Loading Ambulance Stretcher

Item# MS3C-100S

This automatic loading ambulance stretcher is light enough for one person to load.

## Features:

- Constructed of high-strength aluminum alloy, light and durable
- Easily cleaned for sterilization
- 2" thick foamed cushion for extra patient comfort
- Adjustable backrest allows for multiple adjustments keeping the wounded comfortable
- Only one person needed to load the stretcher into the emergency vehicle
- Bumper rails attached to the stretcher legs protect the frame while loading
- Siderails drop to allow for easy patient loading and unloading
- Mechanical foldaway legs are controlled by the handles on either end of the stretcher
- Legs automatically fold under the stretcher when loading into the vehicle
- Locking device included to lock the stretcher in place during transport
- Durable powder-coated frame
- Two Dual locking rear wheels with sealed bearings stabilize patient loading and unloading
- Rotating casters increase the range of motion and sharper turn radius for use in confined spaces
- Light-weight, aluminum construction
- Patient restraints included
- 350 lb load capacity
- 37in maximum load height
- FDA and CE Certified

## Specifications:

- 75 x 22 x 36 highest dimensions
- 75 x 22 x 11 lowest dimensions
- 65-degree maximum backrest angle adjustment

Manufacturer's Warranty: 1 year



Toll Free 1.877.706.4480 | Local/ International 713.706.4480 | Fax 713.706.4161  
[www.MSEC.co](http://www.MSEC.co)