

MS3C-1000PRO

Heavy Duty Folding Lightweight Stretcher-Chair

An affordable, lightweight, durable stretcher that's easy to operate. This stretcher can be folded into a transport chair permitting access to all those confined spaces such as corridors with narrow angles or elevators.

Technical Features:

- Easy to disinfect and clean
- 4-point safety harness included
- I.V. Pole included
- Stretcher can be adjusted to 5 pre-set positions
- Mattress is made of 2in waterproof foamed cushion
- 2 locking dropdown siderails
- Height of the backrest is adjusted with a pneumatic spring
- Made of high-strength, lightweight aluminum alloy material
- Ideal solution for patients who cannot be subjected to long waiting time in a flat position
- Only one person required to operate the stretcher when loading it into an emergency vehicle
- Comes standard with mechanical folding legs that are controlled by handles on both sides and automatically lock in place when the undercarriage is fully unfolded
- Limited Manufacturer Warranty: 1 Year



YouTube Demonstration:

MS3C-1000Pro Heavy Duty Folding Lightweight Stretcher-Chair

MS3C-1000Pro Seat Belt Installation



INCLUDED ACCESSORIES



2in Waterproof
PVC Foam
Mattress



(3) Vest Type
Belts: Chest,
Hip & Knee



I.V.
Pole

Model	Upright Dimensions (L x W x H)	Collapsed Dimensions (L x W x H)	Full Upright Chair Dimensions (L x W x H)	Load Bearing Weight	Maximum Height
					Vehicle Loading
MS3C-1000PRO	77in x 22in x 33.5in	77in x 22in x 34in	77in x 22in x 18in	550lbs	26.5in