

## MS3C-250 PRO Commercial Multi Height Aluminum Alloy EMS Stretcher

### Standard Features:

- Multifunctional automatic loading stretcher has 5 preset heights with conveniently placed releases
- The bed surface is made of 4" contoured foamed cushion and the adjustable back allowing for comfortable patient transport
- 4 Point Chest Harness with belt and (2) Patient Safety belts included
- Adjustable footrest
- Positions include Trendelenburg / Reverse Trendelenburg
- Removable top works as a backboard
- (4) 8in Heavy duty, (2) locking swivel to provide controlled steering and maneuvering
- Locking sidebars collapse to allow for easy patient loading and unloading
- Vehicle locking device included to lock unit in place when transporting patient
- One operator for loading and unloading into transportation vehicle
- Easily cleaned and sterilized
- Adjustable IV Pole
- O2 holder
- Vehicle Loading Height 25in (Distance from the ground to the loading platform ie. bumper)



### Positions (top of pad):

- Position 1: 77 x 22 x 42
- Position 2: 77 x 22 x 36
- Position 3: 77 x 22 x 30
- Position 4: 77 x 22 x 25
- Position 5: 77 x 22 x 20

### Specifications:

- Weight: 145 lbs
- Load capacity: 550 lbs

**Manufacturer warranty: 1 year**



Toll Free 1.877.706.4480 | Local/International 713.706.4480 | Fax 713.706.4161

[www.MS3C.co](http://www.MS3C.co)