



toll free: 1-877-706-4480
international/local: 713-706-4480
fax: 713-706-4161
www.MS3C.co

EMS Lite Evacuation Stair Chair Item # MS3C-330ML

Features:

- Extendable handle allows the operator to have a forward facing, upright posture for balance
- Telescoping foot end lift handles also extend, allowing two people to lift over obstacles or up stairs.
- Easy set-up and instruction guide on seatback
- Heavy duty track allows easily controlled stair descent
- Telescoping head end and foot end lift handles for ergonomic lifting
- Molded hand grips fit hands for better ergonomics
- Positive action locking mechanism for passenger security
- Dual wheel locks for more secure transfers
- Oversized rear wheels with sealed bearings for smoother transports
- Compact storage size conserves space in corridors and closets
- Lightweight, durable, aluminum construction for low effort transports
- Grease-free maintenance for ease of care
- Passenger safety restraints
- Head support allows evacuation of passengers with compromised head control
- Easy to clean two-piece molded ABS seat for passenger comfort
- Available in red or yellow

Specifications:

- Height 45"
- Width 20.5"
- Step height: 4.5"
- Floor clearance: 2.25"
- Weight: 30 pounds
- Max. weight capacity: 500 pounds

